

Sheet No.

**GT200-PE021E Oil**

## Determination of acid number in gasoline engine oil 1/6

\*This application sheet is provided as reference, and does not assure the measurement results. Please consider analysis environment, external factors and sample nature for optimal conditions before the measurement.

### Outline

Acid number in gasoline engine oil is determined with titration by potassium hydroxide in 2-propanol titrant after dissolving new or used oil in titration solvent contains toluene, 2-propanol and small amount of water. The titration result is used as reference of oxidation and deterioration state of the oils for example.

Titration Type : Non-aqueous Neutralization, Titration mode: INF, Detection: pH/mV  
 ◆Reference : **ASTM D664-07** Standard Test Method for Acid Number of Petroleum Products by Potentiometric Titration

### Apparatus

Automatic titrator : GT-200  
 Electrodes : Reference Electrode sleeve type, Glass electrode  
 Reference electrode solution : 3mol/L, Lithium chloride in ethanol  
 Buret size : 10ml

### Reagents

[ Titrant ]

■0.1mol/L- Potassium hydroxide in 2-propanol for testing neutralization number in oil

[ Prepared reagents ]

■Titration solvent : mixed 500ml of toluene, 495ml of 2-propanol and 5ml of pure water

■3mol/L of Lithium chloride in ethanol : Dissolve 12.7g of lithium chloride, special grade reagent, in ethanol, special grade reagent, and dilute the solution to 100ml by the ethanol.

### Analytical Procedure

[ Blank measurement ]

- (1) Add 125ml of the titration solvent into a 200ml beaker by a measuring cylinder.
- (2) Titrate with 0.1mol/L . potassium hydroxide in 2-propanol titrant

[ Sample measurement ]

- (1) Add proper size of sample decided by the method depending on the acid number of the sample into a 200ml beaker. 5g +/-0.5g in this sample.
- (2) Add 125ml of the titration solvent into the above mentioned 200ml beaker by a measuring cylinder.
- (3) Titrate with 0.1mol/L . potassium hydroxide in 2-propanol titrant

Sheet No.

**GT200-PE021E Determination of acid number in gasoline engine oil** ————— **2/6**

[ Calculation ]

**Acid number ( mgKOH/g ) = ( A1 - BL ) × M × E × f × FW / S × R**
**(Used prefixed formula on GT-200)**

A1 : Titration volume of 0.1mol/L- potassium hydroxide in 2-propanol titrant for sample measurement (ml)

BL : Titration volume of 0.1mol/L- potassium hydroxide in 2-propanol titrant for blank measurement (ml)

M : Molarity of 0.1mol/L- potassium hydroxide in 2-propanol titrant (0.1)

E : Equivalent number of 0.1mol/L- potassium hydroxide in 2-propanol titrant (1)

F : Factor of 0.1mol/L- potassium hydroxide in 2-propanol titrant

FW : Formula weight of potassium hydroxide (56.1)

S : Sample size(g)

R : Dilution rate (1)

### Other Requirements

- When measuring samples with pH detection, calibrate the apparatus by three standards, pH 7, 4 and 11 before measurement. Select “Sleeve type liquid: 3.3M KCL (GTRS10B)” and “Three point calibration (Input pH)” on the “pH Calibration” of GT-200.
- For using 10ml Buret, set the volume by “Setting” on the Automatic Buret’s software.
- After a measurement, wash the electrodes by the titration solvent and immerse them in pure water for 5min. as conditioning.
- Confirm reagent labels and safety data sheets for safety
- Wear protective equipment (eye protector, gloves and others.) when handling reagents.

### Measurement Results

Detection : mV

	Sample size(g)	Titrant (ml)	Results(mg KOH/g)
1	5.0218	2.5412	2.71
2	5.0136	2.5079	2.67
3	5.0592	2.5590	2.71

N            3  
 Average    2.69  
 SD           0.019  
 RSD(%)    0.70

Detection : pH

	Sample size(g)	Titrant (ml)	Results(mg KOH/g)
1	5.0617	2.5524	2.69
2	5.0292	2.5016	2.65
3	5.0388	2.5101	2.65

N            3  
 Average    2.66  
 SD           0.022  
 RSD(%)    0.81

Acid number in gasoline engine oil (0W-20) is measured by GT-200.

Average of three measurements is around 2.7mgKOH / g. The results are repeatable on both mV and pH detections.

Sheet No.


**GT200-PE021E** Determination of acid number in gasoline engine oil 3/6

ID No. : 1 GT No.1

User : GT-200

Measurement : 2014/08/19 12:00  
 Sample Name : BLANK

Type : Sample Titr  
 Sample Size : 125 [ml]



C1 : 0.1193 [ml]

A1 : 0.1193 [ml] -168 [mV]

P-initial : -15 [mV]  
 Start : 0 [ml] -15 [mV]  
 End : 0.5 [ml] -495 [mV] Measuring Time : 5'16"

File No. : 14 OIL / Acid Number  
 Titr File No. : 39 Acid Number / Blank  
 Mode : INF End1, End1 Width : -500 [mV] ± 1000 [mV]  
 Detect : mV1  
 BRT No. : 1  
 Reagent : 13  
 WTint : 30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 10 [mV]  
 dT : 10 [sec]  
 DL : 800 [mV/ml]  
 DetCnt : 3  
 Vmax : 20 [ml]  
 Vover : 0.2 [ml] C1 : A1 [ml]  
 Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
 F : 1  
 Buret Injection Speed : 500 [ul/sec]

Sheet No.


**GT200-PE021E** Determination of acid number in gasoline engine oil 4/6

ID No. : 5 GT No.1

User : GT-200

Measurement : 2014/08/19 15:25  
Sample name : Engine Oil

Type : Sample Titr  
Sample Size : 5.0218 [g]



C1 : 2.71 [mgKOH/g]

A1 : 2.5412 [ml] -330 [mV]

P-initial : 20 [mV]  
Start : 0 [ml] 20 [mV]  
End : 3.738 [ml] -479 [mV] Measuring Time : 13'31"

File No. : 14 OIL / Acid Number  
Titr File No. : 6 OIL / Acid Number  
Mode : INF End1, End1 Width : -300 [mV] ± 500 [mV]  
Detect : mV1  
BRT No. : 1  
Reagent : 13  
WTint : 30 [sec]  
Vup : 100 [μl]  
Vlow : 50 [μl]  
dE : 10 [mV]  
dT : 10 [sec]  
DL : 200 [mV/ml]  
DetCnt : 20  
Vmax : 20 [ml] C1 : (A1-BL)\*M\*E\*f\*FW/S\*R  
Vover : 0.2 [ml] [mgKOH/g]

Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
F : 1 BL : 0.1193 [ml]  
FW : 56.1 R : 1

Buret Injection Speed : 500 [ul/sec]

Sheet No.


**GT200-PE021E** Determination of acid number in gasoline engine oil 5/6

ID No. : 1GT No.1

User : GT-200

Measurement : 2014/08/20 10:10  
 Sample name : Blank

Type : Sample Titr  
 Sample size : 125 [ml]



C1 : 0.1261 [ml]

A1 : 0.1261 [ml] 11.229 [pH]

P-initial : 6.939 [pH]  
 Start : 0[ml] 6.939 [pH]  
 End : 1.75 [ml] 14 [pH]

Measuring Time : 6'15"

FileNo. : 14OIL / Acid Number  
 Titr FileNo. : 40 AcidNumber/BlankpH  
 Mode : INF End1, End1 Width : 11 [pH] ± 2 [pH]  
 Detect : pH  
 BRT No. : 1  
 Reagent : 13WTint:30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 0.2 [pH]  
 dT : 10 [sec]  
 DL : 20 [pH/ml]  
 DetCnt : 3  
 Vmax : 20 [ml]  
 Vover : 0.2 [ml] C1 : A1

[ml]

Reag : 0.1M-KOH/IPA E : 1 M : 0.1 [Mol/l]  
 F : 1

BuretInjectionSpeed: 500 [ul/sec]

Sheet No.


**GT200-PE021E** Determination of acid number in gasoline engine oil 6/6

ID No. : 2 GT

No.1 User : GT-200

Measurement : 2014/08/20 10:47  
 Sample name : Engine oil

Type : Sample Titr  
 Sample Size : 5.0617 [g]



C1 : 2.69 [mgKOH/g]

A1 : 2.5524 [ml] 12.521[pH]

P-initial : 6.352 [pH]  
 Start : 0 [ml]  
 End : 3.578 [ml]

6.352 [pH]  
 14 [pH]

Measuring Time : 11'30"

File No. : 14 OIL / Acid Number  
 Titr File No. : 36 Acid Number pH  
 Mode : INF  
 Detect : pH  
 BRT No. : 1  
 Reagent : 13  
 WTint : 30 [sec]  
 Vup : 100 [μl]  
 Vlow : 50 [μl]  
 dE : 0.2 [pH]  
 dT : 10 [sec]  
 DL : 5 [pH/ml]  
 DetCnt : 10  
 Vmax : 20 [ml]  
 Vover : 0.2 [ml]

End1, End1 Width : 11 [pH] ± 2 [pH]

C1 : (A1-BL)\*M\*E\*f\*FW/S\*R

[mgKOH/g]

Reag : 0.1M-KOH/IPA  
 F : 1  
 FW : 56.1

E : 1  
 BL : 0.1261 [ml]  
 R : 1  
 M : 0.1 [Mol/l]

Buret Injection Speed : 500 [ul/sec]