

Sheet No.

**GT200-FO001** Vinegar

## Determination of Acidity in fermented vinegar ( grain vinegar ) \_\_\_\_\_ 1/3

\*This application sheet is provided as reference, and does not assure the measurement results. Please consider analysis environment, external factors and sample nature for optimal conditions before the measurement.

### Outline

Determination of Acidity in fermented vinegar is identified in the Japan Agricultural Standards for fermented vinegar. It provides that the acidity of grain vinegar must be 4.2% or more.

Titration Type : Neutralizing  
◆Reference : JAS for fermented vinegar  
Acidity determination : Automatic titration(method for potentiometric titrator)

### Apparatus

Automatic titrator : GT-200  
Electrodes : Double junction type reference electrode, Glass electrode  
Reference electrode solution : Inner : 1 mol / L - potassium chloride in water  
Outer : 1 mol / L - potassium nitrate in water

### Reagents

[ Titration Solution ] ■0.5mol / L - Sodium hydroxide in water ( Volumetric analysis grade )

### Analytical Procedure

[ Blank measurement ]

- (1) Add 100ml pure water into a 200ml beaker by measuring cylinder.
- (2) Titrate with 0.5 mol / L - Sodium hydroxide solution. ( MODE : SET-P, END 1 : 8.2 pH )

[ Sample measurement ]

- (1) Add sample into a 200ml beaker by volumetric pipette. Adjust sample volume so that the titrant consumption will be 10 to 20 ml. (Sample volume of this application sheet is 10 ml.)
- (2) Add 100 ml pure water into a beaker by measuring cylinder.
- (3) Titrate with 0.5mol / L - Sodium hydroxide solution. ( MODE : SET-P, END 1 : 8.2 pH )

Sheet No.

**GT200-FO001 Determination of Acidity in fermented vinegar ( grain vinegar )** ——— 2/3

[ Calculation ]

$$\text{Acidity (\%)} = 0.03 \times (\text{A1-BL}) \times f / \text{S} \times 100$$

0.03 : Weight of acetic acid equivalent to 1ml of 0.5 mol / L - Sodium hydroxide solution (g)

A1 : Titration volume of 0.5 mol / L - Sodium hydroxide solution for Sample measurement (ml)

BL : Titration volume of 0.5 mol / L - Sodium hydroxide solution for blank measurement (ml)

F : Factor of 0.5 mol / L - Sodium hydroxide solution

S : Sample Volume (ml)

**Other Requirements**

- pH calibration with pH standard solution is required before measurement.
- Confirm reagent labels and safety data sheets for safety.
- Wear protective equipment (eye protector, gloves and others.) when handling reagents.

**Measurement Results**

	Sample Volume (ml)	Titration volume (ml)	Measurement value (%)
1	10	14.1230	4.2
2		14.1189	4.2
3		14.1243	4.2

N                    3  
 Average           4.2  
 SD                   0.0008  
 RSD (%)           0.0200

The result shows average of three times measurement is 4.2% and RSD is 0.02%. GT-200 can measure determination of acidity in fermented vinegar ( grain vinegar ) with good repeatability.

Sheet No.


**GT200-FO001** Determination of Acidity in fermented vinegar ( grain vinegar ) 3/3

ID No. : 3 GT No.1

User : GT-200

Measurement : 2013/03/11 10:39  
 Sample Name : Grain vinegar

Type : Sample Titr  
 Sample Size : 10 [ml]



C1 : 4.24 [%]

A1 : 14.1230 [ml] 8.2 [pH]

P-Initial : 3.177 [pH]  
 Start : 0 [ml] 3.177 [pH]  
 End : 14.704 [ml] 11.65 [pH] Measuring Time : 5'18"

File No. : 0 Quick Mode  
 Titr File No. : 16 Acidity determination of fermented vinegar (grain vinegar)  
 Mode : SET-P End1 : 8.2[pH]  
 Detect : pH  
 BRT No. : 1  
 Reagent : 12  
 WTint : 0 [sec]  
 Vup : 400 [μl]  
 Vlow : 10 [μl]  
 dE : 0.1 [pH]  
 dT : 3 [sec]  
 Vmax : 25 [ml]  
 Vover : 0.5 [ml]

$$C1 : 0.03 \cdot (A1 - BL) \cdot f / S \cdot 100 \quad [\%]$$

Reag : NaOH E : 1 M : 0.5 [Mol/l]  
 F : 1.002 BL : 0.0049 [ml]

Buret Injection Speed : 500 [ul/sec]