

Sheet No.

**AQF PE 006E** Oil

## Determination of fluorine in lubricating oil ————— 1/2

Instruments : AQF-100

Method : Combustion-ion chromatography

Related standard :

Concentrations of fluorine, chlorine, bromine, iodine, and sulfur can be determined and accurately by using a combustion ion chromatography (CIC) system combining an Automatic Quick Furnace Model AQF-100 which safely combusts samples with an ion chromatograph.

|                       |   |
|-----------------------|---|
| Sample name           | Lubricating oil   |
| Sample status         |   |
| Measuring items       | Fluorine (F)  |
| Measurement principle | <p>Sample is thermally decomposed in argon (Ar) atmosphere, then combusted in oxygen (O<sub>2</sub>) atmosphere. Halogens in the sample are converted to hydrogen halide and halogen gas and sulfur turns into sulfur oxide. These components are collected into absorbing solution and converted to halide ion and sulfate ion. The resulting solution is analyzed by injecting into an ion chromatograph (IC).</p> <p><b>Analyzing flow</b><br/>         [Sample weighing]→[Combustion]→[Collection of combustion gas]→[IC analysis]</p>  |
| Parameters            | <p><b>1. AQF-100</b></p> <p>Sample size : 50mg<br/>         Sample boat : Quartz sample boat, TX2SBT<br/>         Additive : Not used<br/>         Pyrolysis tube : Quartz tube filled with quartz wool<br/>         Absorbent : Hydrogen peroxide / water<br/>         Mode :</p> <p>Heater Temp. Inlet : 800degC<br/>         Outlet : 1000degC<br/>         Gas flow Ar : 200 ml/min<br/>         O<sub>2</sub> : 400 ml/min</p> <p>GA-100 Absorbent volume : 5ml<br/>         Sampling loop : 20 μl<br/>         Absorption tube : For 10 ml<br/>         Water supply : 2<br/>         Ar flow for water supply : 150 ml/min</p> |

Sheet No.

**AQF PE 006E Determination of fluorine in lubricating oil**


2/2

## 2. Ion chromatograph

Ion chromatograph : DIONEX DX-120  
 Column : DIONEX Ion Pack AG12A / Ion Pack AS12A  
 Eluent : 2.7mM Na<sub>2</sub>CO<sub>3</sub> / 0.3mM NaHCO<sub>3</sub>  
 Eluent flow : 1.50ml / min  
 Detector : Conductivity  
 Suppressor : SRS  
 Measuring time : 15min  
 Sampling loop : 100 μl using GA-100 sampling loop  
 Calibration : F Cl Br S :0.1ppm to 5.0ppm

Results

### Chromatogram



### Results

| Sample   | Result (ppm) | Average (ppm) |
|----------|--------------|---------------|
| Sample A | 2.5 , 2.7    | 2.6           |
| Sample B | 10.5 , 10.3  | 10.4          |

Remarks

\*Handling of reagents: Confirm labels and safety data sheets of reagents and handle them with enough care.  
 \*Automation is possible by using an Automatic Sample Changer, ASC-120S.  
 \*When ASC-120S is used, the boat to be used will be a ceramic boat, TX3SCX.

\*This application sheet is provided as reference, and does not assure the measurement results. Please consider analysis environment, external factors and sample nature for optimal conditions before the measurement.

AQF100\_03\_004E