

## NIST No. 1083 (Wear-Metals in Lubricating Oil)

Sheet No.: Measurement model: NSX2100H-CH-001E TS-2100H System ABC-210/HF-210/SD-210

Relevant standard:

Petroleum chemistry ASTM D5453 JIS K 2541

Detection method: Ultraviolet Fluorescence method

The sulfur in lubrication oil is added as extreme pressure agent. It needs to be managed at low concentration from a standpoint of environmental load. The sulfur analysis device (**TS-2100H**) of Mitsubishi Chemical Analytech Co., Ltd. can analyze the sulfur in lubrication oil quickly with accuracy.

Sample name	Wear-Metals in Lubricating Oil								
Analytical item	Quantitative analysis of sulfur in combustion method								
Standard	<u>ASTM-D5453</u> : standard testing method for measuring sulfur contained in carbon hydride and fuel using an ultraviolet fluorescence detector <u>JIS K 2541</u> : raw petroleum and petroleum product – sulfur content testing method – ultraviolet fluorescence method								
Analytical principle Result of sulfur	Ultraviolet fluorescence method: Sample is burned in argon / oxygen stream and the generated sulfur dioxide is introduced to the cell of ultraviolet irradiation. The fluorescence intensity generated by ultraviolet irradiation is measured and the amount of sulfur is calculated based on the standard curve that has been created using the standard sulfur sample. Organic-S + O <sub>2</sub> $\rightarrow$ SO <sub>2</sub> + CO <sub>2</sub> (combustion) SO <sub>2</sub> + hv $\rightarrow$ SO <sub>2</sub> + hv <sub>2</sub> (ultraviolet fluorescence)								
analysis Horizontal combustion method	Comple nome	<b>TS-2100H</b> analysis value (S μg/g)							
	Sample name	1	2	3	Average	RSD (%)			
	Wear-Metals in Lubricating Oil	1,017	1,016	1,016 1,014		0.15			
Required analysis time	Pretreatment () minutes, Measurement ( 9) minutes Total ( 9) minutes/ (1) measurement								
Horizontal type									

\*This sheet is provided as a reference and does not guarantee analytical values. Optimal conditions may vary depending on external factors, such as the analysis environment, and the nature of the sample.



	Measurement conditionGas flow rateReaction tube double tube for ABCAr										
O an dition of	[ABC program]		2nd		3rd		End	Cool	Delay	1	
Conditions of sulfur analysis	Sample name	Pos (mm)	Time (sec)	Pos (mm)	Time (sec)	Pos (mm)	Time (sec)	Time (sec)	Time (sec)	Time (sec)	
Horizontal combustion	Lubrication oil	75	10	95	120	105	20	100	60	180	
method	hod Boat Speed: 20mm/sec Ar Time: 10sec O <sub>2</sub> Time (sec): 600sec								5		
	<ul> <li>Standard sample for standard curve: S_Dibutyl disulfide / toluene 0, 5, 50, 250µg/mL × 40µL</li> <li>Amount of introduced sample: 20µL</li> <li>For the sample for measurement, the sample diluted to 10% with toluene used. The obtained measurement value multiplied by dilution rate was set to sulfur quantitative value.</li> </ul>										