## Nittoseiko Analytech



## NIST No. 1616b (Sulfur in Kerosene (Low Level))

Sheet No.: NSX2100H-PE-003E Petroleum chemistry

Measurement model: TS-2100H System Relevant standard: ASTM D5453 ABC-210/HF-210/SD-210 JIS K 2541

Detection method: Ultraviolet Fluorescence method

The sulfur content in heating oil becomes sulfur oxide (SOx) during combustion and directly affects the atmospheric air pollution as substance of concern. For this reason, it needs to be low concentration and the quantitation method complying it is required. The sulfur analysis device (**TS-2100H**) of Mitsubishi Chemical Analytech Co., Ltd. can analyze the sulfur in heating oil quickly with accuracy.

Sample name	Sulfur in Kerosene (Low Level)									
Analytical item	Quantitative analysis of sulfur in combustion method									
Standard	ASTM-D5453: standard testing method for measuring sulfur contained in carbon hydride and fuel using an ultraviolet fluorescence detector  JIS K 2541: raw petroleum and petroleum product – sulfur content testing method – ultraviolet fluorescence method									
Analytical principle  Result of sulfur	Ultraviolet fluorescence method: Sample is burned in argon / oxygen stream and the generated sulfur dioxide is introduced to the cell of ultraviolet irradiation. The fluorescence intensity generated by ultraviolet irradiation is measured and the amount of sulfur is calculated based on the standard curve that has been created using the standard sulfur sample. $ Organic-S + O_2 \rightarrow SO_2 + CO_2  (combustion) \\ SO_2 + hv \rightarrow SO_2 + hv_2  (ultraviolet fluorescence) $									
analysis Horizontal combustion method		TS-2100H analysis value (S mg/kg)								
	Sample name	1	2	3	Average	RSD (%)				
	Sulfur in Kerosene (Low Level)	8.28	8.21	8.19	8.23	0.55				
Danisad										
Required analysis time	Pretreatment () minutes, Measurement ( 9) minutes Total ( 9) minutes/ (1) measurement									
Horizontal type	istal (5) illinoista (7) illoudatoriista									

https://www.n-analytech.co.jp/

<sup>\*</sup>This sheet is provided as a reference and does not guarantee analytical values. Optimal conditions may vary depending on external factors, such as the analysis environment, and the nature of the sample.

## Nittoseiko Analytech



Gas flow rate

Ar...... 300mL/min

O<sub>2</sub> ...... 300mL/min

O<sub>2</sub> Time (sec): 600sec

## **Measurement condition**

Reaction tube ... double tube for ABC Temperature of electric furnace

Inlet Temp 800°C Outlet Temp 1,000°C

Boat Speed: 20mm/sec

PMT Range High (for Low concentration)

Conditions of

sulfur analysis

Horizontal combustion method

[ABC program]

ABC program												
Sample	1st		2nd		3rd		End Time	Cool Time	Delay Time			
name	Pos (mm)	Time (sec)	Pos (mm)	Time (sec)	Pos (mm)	Time (sec)	(sec)	(sec)	(sec)			
Lubrication oil	105	20	125	80	140	20	100	60	200			

Ar Time: 30sec

Standard sample for standard curve: S\_Dibutyl disulfide / toluene 0, 1, 10,  $50\mu g/mL \times 20\mu L$ 

Amount of introduced sample: 20µL

- The sample was not diluted, but was introduced directly.