

## NIST No. 1083 (Wear-Metals in Lubricating Oil)

Sheet No.: Measurement model: NSX2100V-CH-002E TS-2100V System ASC-250L/VF-210/SD-210

Relevant standard:

Petroleum chemistry ASTM D5453 JIS K 2541

Detection method: Ultraviolet Fluorescence method

The sulfur in lubrication oil is added as extreme pressure agent. It needs to be managed at low concentration from a standpoint of environmental load. The sulfur analysis device (**TS-2100V**) of Mitsubishi Chemical Analytech Co., Ltd. can analyze the sulfur in lubrication oil quickly with accuracy.

Sample name	Wear-Metals in Lubricating Oil					
Analytical item	Quantitative analysis of sulfur in combustion method					
Standard	ASTM-D5453: standard testing method for measuring sulfur contained in carbon hydride and fuel using an ultraviolet fluorescence detector JIS K 2541: raw petroleum and petroleum product – sulfur content testing method – ultraviolet fluorescence method					
Analytical principle Result of sulfur analysis Vertical combustion method	Ultraviolet fluorescence method: Sample is burned in argon / oxygen stream and the generated sulfur dioxide is introduced to the cell of ultraviolet irradiation. The fluorescence intensity generated by ultraviolet irradiation is measured and the amount of sulfur is calculated based on the standard curve that has been created using the standard sulfur sample. Organic-S + O <sub>2</sub> $\rightarrow$ SO <sub>2</sub> + CO <sub>2</sub> (combustion) SO <sub>2</sub> + hv $\rightarrow$ SO <sub>2</sub> + hv <sub>2</sub> (ultraviolet fluorescence)					
	Sample name	<b>TS-2100V</b> analysis value (S µg/g)				
		1	2	3	Average	RSD (%)
	Wear-Metals in Lubricating Oil	1,016	1,022	1,033	1,024	0.83
Deguired						
Required analysis time	Pretreatment () minutes, Measurement ( 4) minutes Total ( 4) minutes/ (1) measurement					
Vertical type						

\*This sheet is provided as a reference and does not guarantee analytical values. Optimal conditions may vary depending on external factors, such as the analysis environment, and the nature of the sample.



	<b>Measurement condition</b> Temperature of electric furnace	Gas flow rate Ar 100mL/min O <sub>2</sub> 500mL/min			
	Inlet Temp 900°C	Ar time: 30sec			
Conditions of sulfur analysis	Outlet Temp 1,000°C PMT Range Low (for High concentration	O <sub>2</sub> time: 120sec			
Vertical combustion method	Standard sample for standard curve: S_Dibutyl disulfide / toluene 0, 1, 10, 50, $100\mu g/mL \times 30\mu L$ Amount of introduced sample: $30\mu L$				
	<ul> <li>For the sample for measurement, the sample diluted to 5% with toluene was used. The obtained measurement value multiplied by dilution rate was set to the sulfur quantitative value.</li> </ul>				