Nittoseiko Analytech

Application Sheet

Sheet No. AQF2100 IM 048 Inorganic Materials

Determination of fluorine in metal fluorides by combustion ion chromatography

Instruments : AQF-2100H system Detection method : Ion chromatography Related standard :

Introduction

Concentrations of fluorine (F), chlorine (Cl), bromine (Br), iodine (I), and sulfur (S) can be determined accurately by combustion ion chromatography (CIC), a system in which samples are safely combusted by Automatic Quick Furnace (model AQF-2100H) and injected online to an ion chromatograph (IC).

Method

Combustion ion chromatography

Sample was pyrolysed in argon (Ar) carrier gas, then pyrolysis gas was combusted in oxygen (O₂) gas. Halogens in the sample were converted into hydrogen halide (HX) and/or halogen gas (X₂), and sulfur in the sample was converted into sulfur oxide (SO_X), which were collected and further converted into halide ion (X⁻) and/or sulfate ion (SO₄²⁻) in the absorbent. This absorbent was injected to IC, and analyzed.

Flowchart


Weighing \Rightarrow combustion \Rightarrow collection \Rightarrow IC


Samples

Metal fluorides

Aluminum fluoride (AIF₃), Yttrium fluoride (YF₃), Magnesium fluoride (MgF₂), Barium fluoride (BaF₂)

Results





*This sheet is provided as a reference and does not guarantee analytical values. Optimal conditions may vary depending on external factors, such as the analysis environment, and the nature of the sample.